

# RTZ-B100 DELUXE SWIVEL SHOWER CHAIR

The Swivel Shower Chair makes showering and cleaning easier, enabling the client to turn around without twisting their body.

## FEATURES

- 360° rotation with 90° lock intervals
- Seat, back and armrests are padded for comfort
- Armrests can be raised to allow easier access from the side
- Height adjustable
- Sturdy lightweight aluminium frame

## SPECIFICATIONS

Seat Height:	45.5 - 55cm
Seat Width:	40cm
Seat Depth:	33cm
Seating Width:	50cm
Seating Depth:	37cm
Backrest Height:	33.5cm
Overall Width:	60cm
Overall Depth:	54cm
Weight Capacity:	150kg
Net Weight:	6.4kg

Assembly required

# RTZ-B101 BARIATRIC ALUMINIUM SHOWER CHAIR



Ideal for people who are unable to stand for long periods of time. Shower chairs enable a person to bathe safely while remaining seated. It gives the user ongoing comfort, stability and greatly reduces the chance of slipping and falling.

## FEATURES

- Lightweight aluminium with powdercoat finish
- Wide 56cm seating space
- Side cross bars to assist in lowering onto and rising from seat
- Height adjustable in 25.4mm (1") increments
- Drainage holes in seat reduce water build up
- Non-slip tips – slip resistant

## SPECIFICATIONS

Seat height range:	41.5 - 56.5cm
Seat width:	44cm
Seating depth:	41cm
Width between arms:	56cm
Armrest height:	29cm
Width at base:	74cm
Depth at base:	51cm
Net weight:	4.5kg
Weight capacity:	250kg

# RTZ-B200 DELUXE FOLDABLE COMMUNE SHOWER CHAIR



A commode and shower chair which can fold neatly for transport or storage where space is limited

## FEATURES

- Seat, back and armrests are padded for comfort
- Armrests can be raised to allow easier access from the side
- Height adjustable
- Sturdy lightweight aluminium frame

## SPECIFICATIONS

Seat Height:	49 - 58.5cm
Seat Width:	41cm
Seat Depth:	40.5cm
Seating Width:	51cm
Seating Depth:	42cm
Backrest Height:	37cm
Armrest Height:	11.5cm
Under Seat Clearance:	37cm - 47cm
Inner Leg Clearance:	42cm
Overall Width	66.5cm
Overall Depth:	55cm
Weight Capacity:	150kg
Folded Dimensions:	60 x 50 x 30cm
Net Weight:	5.8kg

Assembly required

# RTZ-B300 HEAVY DUTY 3 in 1 SHOWER CHAIR COMMODE



A versatile heavy duty commode which can also be used as an over toilet aid and shower chair.

## FEATURES

- Seat, back and armrests are padded for comfort
- Height adjustable
- Sturdy lightweight aluminium frame

## SPECIFICATIONS

Seat Height:	49 - 59cm
Seat Width:	46cm
Seat Depth:	41cm
Seating Width:	57cm
Seating Depth:	45cm
Backrest Height:	32.5cm
Armrest Height:	18cm
Under Seat Clearance:	39 -49cm
Inner Leg Clearance:	50cm
Overall Width:	71cm
Overall Depth:	55cm
Weight Capacity:	250kg
Net Weight:	5.5kg

Assembly required

# RTZ-B400 PADDED ROTATING SLIDING BATH SEAT



A Transfer Bath Bench with a rotating and sliding seat makes transfers to and from a bath tub easy.

## FEATURES

- 360° rotation with 90° lock intervals
- Seat slides along rails to position over tub
- Padded seat, back and armrests for comfort
- Support base can be adjusted to suit most standard bath widths.
- Sturdy lightweight aluminium frame

## SPECIFICATIONS

Seat Height	14cm (from support base)
Seat Width:	40cm
Seat Depth:	32.5cm
Seating Width:	54cm
Seating Depth:	36cm
Backrest Height:	29cm
Armrest Height:	15cm
Fits Tub Width:	56 - 71cm (2.5cm increments)
Overall Width:	70.5 - 85.5cm
Overall Depth:	55.5cm
Weight Capacity:	150kg
Net Weight:	8.25kg

Assembly required

# RTZ-E100 ARM & LEG EXERCISER

Ideal for users with limited movement and/or requiring arm and leg rehabilitation. It enables low impact exercise to help tone muscles and increase blood circulation.

## FEATURES

- Adjustable resistance
- Angle of handles can be adjusted to suit reach
- Steel tube construction

## SPECIFICATIONS

Length:	54cm
Width:	42cm
Height:	98cm
Net weight:	6.75kg
Max. loading weight:	40kg

Assembly required



© 2022 Pharmaquip Healthcare

# RTZ-T100 BARIATRIC OVER TOILET AID

Over Toilet Frames are designed to increase the height of the toilet seat and includes arm supports to assist the user in transferring on and off the toilet safely.

## FEATURES

- White powdercoated lightweight aluminium frame
- Height adjustable in 25mm (1") increments
- Wide 56cm seating width for a comfortable seating space
- Supplied with easy to clean removable seat, lid and splash guard
- Seat and lid made from high density plastic to minimise flex
- Side bars to assist in lowering onto and rising from seat
- Non-slip tips – slip resistant
- Conditional warranty applies

## SPECIFICATIONS

Seat height range:	42-57cm
Width between arms:	56cm
Armrest height:	26.5cm
Seating width:	56cm
Seating depth:	42cm
Width at base:	73cm
Depth at base:	51cm
Net weight:	4.8kg
Weight Capacity:	250kg

